



### 1816. Rak koppling

Art.nr	För slang i.d	Hål
1816.902.5.2	4,5	2
1816.902.2	6,0	2
1816.903.7.2	6,8	3
1816.903.2	8,0	3
1816.904.9.2	9,0	4
1816.904.2	10,0	4
1816.905.2	11,0	5



### 1816. L-Koppling

Art.nr	För slang i.d	Hål
1816.1002.5.2	4,5	2
1816.1002.2	6,0	2
1816.1003.7.2	6,8	3
1816.1003.2	8,0	3
1816.1004.9.2	9,0	4
1816.1004.2	10,0	4
1816.1005.2	11,0	5



### 1816. NS-koppling

Art.nr	För slang i.d	NS
1816.2201.2	8	14
1816.1005.2	8	19

Kopplingarna är i ren PTFE och kan steriliseras vid 135°C



### 1816. Backventil

Art.nr	För slang i.d	Hål
1816.2402.5.2	4,5	2
1816.2403.2	6	2
1816.2404.9.2	9	4
1816.2404.2	10	4

• Mått i mm •