


**1816. Backventil**

Art.nr	För slang y.d	Hål
1816.2302.6.2	6	2
1816.2304.2	8	4
1816.2306.2	12	6
1816.2308.2	14	8

Klarar 4 bar vid 100°C och 1mm vacuum
Kan steriliseras i 135°C

**1816. Rak ventil**

Art.nr	För slang y.d	Hål
1816.1402.6.2	6	2
1816.1403.2	8	3
1816.1404.2	8	4
1816.1406.12.2	12	6
1816.1408.2	14	8
1816.1410.2	14	10

**1816. L-Koppling**


Art.nr	För slang y.d	Hål
1816.1502.6.2	6	2
1816.1503.2	8	3
1816.1504.2	8	4
1816.1506.12.2	12	6
1816.1508.2	14	8
1816.1510.2	14	10

**1816. Reduktionskoppling**

Art.nr	För slang y.d	Hål
1816.1801.2	6	8
1816.1802.2	8	14

**1816. NS-koppling**

Art.nr	För slang y.d	NS
1816.1901.2	8	8
1816.1902.2	8	14

**1816. T-Koppling**

Art.nr	För slang y.d	Hål
1816.1602.6.2	6	2
1816.1603.2	8	3
1816.1604.2	8	4
1816.1606.12.2	12	6
1816.1608.2	14	8
1816.1610.2	14	10