



207



100 ml 150 ml 180 ml 200 ml



250 ml 380 ml 400 ml 500 ml låg



500 ml hög 800 ml

207. Burkar av transparent polypropylen med snäpplock av polyeten

Art.nr.	Volym	Öd	Ud	Höjd	Antal/fp
207.0100PPA	100	61	53	41	2160
207.0150PPA	150	84	76	35	1140
207.0180PPA	180	84	73	42	1020
207.0200PPA	200	84	73	48	1120
207.0250PPA	250	84	68	60	1000
207.0380PPA	380	84	63	95	900

Lock beställes separat

207. Burkar av transparent polypropylen med snäpplock av polyeten

Art.nr	Volym	Öd	Ud	Höjd	Antal/fp
207.0030PS	30	41	35	28	3670
207.0035PS	35	65	64	13	1100
207.0040PS	40	41	34	37	2600
207.0050PS	50	41	32	48	2520
207.0060PS	60	65	65	22	750
207.0100PS	100	65	58	38	1150
207.0170PS	170	65	53	77	750
207.0170PSG	170	65	53	70	775*
207.0205PS	205	65	52	77	100
207.0250PS	250	65	62	84	190
207.0330PS	330	98	85	51	600

*Material slagttålig styren

Lock beställes separat!

207



30 ml 35 ml 40 ml 50 ml



60 ml 100 ml 170 ml 170 ml*



205 ml 250 ml 280 ml 330 ml

• Mått i mm • Volym i ml