

Elektroniska omrörare



more colours on page 20

Stirrers	RW 11 basic	RW 14 basic	RW 16 basic	EUROSTAR digital
Stirring quantity (H ₂ O)	2 l	8 l	10 l	20 l
Viscosity range	100 mPas	10.000 mPas	10.000 mPas	10.000 mPas
Motor rating input / output	8 / 1 W	28,5 / 17 W	75 / 55 W	75 / 55 W
Output at stirring shaft	1 W	17 W	53 W	53 W
Max. ON-time	100 %	100 %	100 %	100 %
Max. torque at chuck (100 rpm)	0,8 Ncm	8 Ncm	40 Ncm	30 Ncm
Speed range I (50 Hz)	0 - 2.000 rpm	100 - 2.000 rpm	40 - 1.200 rpm	50 - 2.000 rpm
Speed range II (50 Hz)	—	—	—	—
Speed display	—	scale	scale (1 - 10)	digital
Chuck range	plug-in coupling 4 mm	0,5 - 8 mm	0,5 - 10 mm	0,5 - 10 mm
Diameter of extension arm	integrated	13 mm / 160 mm	13 mm / 160 mm	13 mm / 160 mm
General information				
Dimensions W x D x H	86 x 175 x 89 mm	70 x 176 x 197 mm	80 x 190 x 222 mm	80 x 190 x 222 mm
Weight	0,39 kg	2 kg	2,8 kg	2,8 kg
Protection class according to DIN EN 60529	IP 42	IP 20	IP 42	IP 42
Interface	—	—	—	—
Torque display	—	—	—	—
Recommended accessories (see page 18)	—	No. 3, 4, 5, 7	No. 3, 4, 5, 7	No. 3, 4, 5, 7
Ident-No.	see page 20	140.3331400	140.2572100	140.2482000