

223

**223. Flaska av aluminium med tätt lock**

Art.nr.	Volym	Diameter	Höjd
223.0060	60	40	88
223.0120	120	57	97
223.0300	300	57	160
223.0600	600	74	192
223.1200	1200	88	256

214

**214. Burk, PE, rund med vid hals, lock av PP**

Art.nr.	Volym	Hals	Diameter	Gradering	Höjd	Antal / fp
214.0050	50	24	40	10	77	400/100
214.0100	100	24	48	20	89	700/100
214.0250	250	35	61	25	126	300/50
214.0500	500	35	75	50	158	160
214.1000	1000	45	93	100	200	80
214.2000	2000	45	115	100	247	40

Levereras med innerlock.

Plomberingsbar.

216

**216. Burk, PE, rektangulär med vid hals, lock av PP**

Art.nr.	Volym	Hals	Basmått	Gradering	Höjd	Antal / fp
216.0025	25	18	32 x 34	5	52	500
216.0050	50	24	38 x 38	10	70	400
216.0100	100	35	42 x 48	20	90	600
216.0250	250	35	57 x 60	25	110	300/50
216.0500	500	45	70 x 80	50	138	150/25
216.1000	1000	58	80 x 103	100	176	80/20
216.2000	2000	58	100 x 134	100	210	40

Levereras med innerlock.

Plomberingsbar.

• Mått i mm • Volym i ml