

337



334



335

**337. Hink PEH**

Art.nr.	Volym	Öd	Ud	Antal/fp	Art.nr lock
337.10	10	270	260	20	337.91
337.12	12	280	310	10	337.92
337.16	16	350	250	10	

12 liters har stålgrepe

334. Emballagehink, av PP

Art.nr.	Volym	Öd	Ud	Höjd	Antal/fp
334.03	3	194	173	135	60
334.05	5	223	200	192	20
334.12	12	297	296	218	15
334.15	15	297	266	292	15
334.20	20	297	255	395	15
339.90					Lock till 3 liter
334.91					Lock till 5 liter
334.92					Lock till 10, 15 och 20 liters hink

335. Hink PEL med plastgrepe

Art.nr.	Volym	Öd	Ud	Höjd	Antal/fp
335.02	2	155	125	170	10
335.03	3	155	135	215	10
335.05	5	200	160	240	10
335.91					Lock till 2 och 3 liter
335.92					Lock till 5 liter

339. Hink PP med pip tål auto-klavering graderad

Art.nr.	Volym	Öd	Höjd	Antal/fp
339.12	12	310	330	1
339.15	15	310	370	1

339

