

**1816. Reduktions Koppling**

| Art.nr | För slang i.d | Hål |
|-------------|---------------|-----|
| 1816.2101.2 | 6 | 8 |
| 1816.2102.2 | 8 | 10 |

**1816. T-Koppling**

| Art.nr | För slang i.d | Hål |
|---------------|---------------|-----|
| 1816.1102.5.2 | 4,5 | 2 |
| 1816.1102.2 | 6,0 | 2 |
| 1816.1103.7.2 | 6,8 | 3 |
| 1816.1103.2 | 8,0 | 3 |
| 1816.1104.9.2 | 9,0 | 4 |
| 1816.1104.2 | 10,0 | 4 |
| 1816.1105.2 | 11,0 | 5 |

**1816. 4-vägs Koppling**

| Art.nr | För slang i.d | Hål |
|---------------|---------------|-----|
| 1816.1112.5.2 | 4,5 | 2 |
| 1816.1112.2 | 6,0 | 2 |
| 1816.1113.7.2 | 6,8 | 3 |
| 1816.1113.2 | 8,0 | 3 |
| 1816.1114.9.2 | 9,0 | 4 |
| 1816.1114.2 | 10,0 | 4 |
| 1816.1115.2 | 11,0 | 5 |

High performance valves

**1816. Rak Koppling med skruv**

| Art.nr | För slang y.d | Hål |
|----------------|---------------|-----|
| 1816.2602.6.2 | 6 | 2 |
| 1816.2603.2 | 8 | 3 |
| 1816.2604.2 | 8 | 4 |
| 1816.2606.12.2 | 12 | 6 |
| 1816.2608.2 | 14 | 8 |
| 1816.2610 | 14 | 10 |

**1816. Rak Koppling**

| Art.nr | För slang i.d | Hål |
|---------------|---------------|-----|
| 1816.2502.5.2 | 4,5 | 2 |
| 1816.2502.2 | 6 | 2 |
| 1816.2503.7.2 | 6,8 | 3 |
| 1816.2503.2 | 8 | 3 |
| 1816.2504.2 | 10 | 4 |
| 1816.2505.2 | 11 | 5 |

Klarar 4 bar vid 100°C och
0,5 mmHg

Alla delar som kommer i
kontakt med slangen är PTFE
och kan steriliseras i 135°C