

# Magnetomrörare med värme



RET basic C safety control

RET basic safety control

RET control-visc C safety control

RET control-visc safety control

2 separate temperature sensors (PT 1000) for heat transfer fluid and medium available (automatic identification)

RET basic safety control	RET basic C safety control	RET control-visc safety control	RET control-visc C safety control
20 l	20 l	20 l	20 l
12 / 5 W	12 / 5 W	12 / 5 W	12 / 5 W
scale	scale	digital	digital
0 - 1.500 rpm	0 - 1.500 rpm	0 - 1.500 rpm	0 - 1.500 rpm
80 x 10 mm	80 x 10 mm	80 x 10 mm	80 x 10 mm
600 W	600 W	600 W	600 W
RT - 340 °C	RT - 340 °C	RT - 340 °C	RT - 340 °C
± 10 K	± 10 K	± 1 K	± 1 K
50 - 380 °C	50 - 380 °C	2 x 50 - 350 °C	2 x 50 - 350 °C
± 1 K	± 1 K	± 0,2 K	± 0,2 K
stainless steel	white coated	stainless steel	white coated
Ø = 135 mm	Ø = 135 mm	Ø = 135 mm	Ø = 135 mm
160 x 280 x 97 mm	160 x 280 x 97 mm	160 x 280 x 97 mm	160 x 280 x 97 mm
2,4 kg	2,4 kg	2,8 kg	2,8 kg
IP 42	IP 42	IP 42	IP 42
—	—	RS 232 / analog	RS 232 / analog
No. 1	No. 1	—	—
140.3188800	140.3197600	140.3364000	140.3364100