


Magnetomrörare utan värme



KMO 2 basic	REO basic C	Midi MR 1 digital	Maxi MR 1 digital
5 l	20 l	50 l	150 l
14 / 4 W	12 / 5 W	70 / 19 W	80 / 35 W
scale	scale	digital	digital
0 - 1.100 rpm	0 - 1.500 rpm	0 - 1.000 rpm	0 - 600 rpm
—	—	—	—
50 x 8 mm	80 x 10 mm	80 x 10 mm	155 x 27 mm
stainless steel	white coated	stainless steel	stainless steel
140 x 120 mm	Ø = 135 mm	350 x 350 mm	500 x 500 mm
140 x 200 x 75 mm	160 x 280 x 97 mm	360 x 430 x 110 mm	505 x 585 x 110 mm
1,4 kg	2,4 kg	9 kg	16 kg
IP 21	IP 42	IP 21	IP 21
140.2812000	140.3261100	140.2621900	140.2621800