

Magnetomrörare utan värme

176

Order order@noaxlab.se



Magnetic stirrers	Mini MR standard	lab disc	COLOR SQUID	BIG SQUID
Stirring quantity (H ₂ O)	0,8 l	0,8 l	0,8 l	1,5 l
Motor rating input / output	4 / 0,5 W	5 / 3 W	2 / 1 W	2 / 1 W
Speed display	—	—	—	—
Speed range	0 - 1.500 rpm	15 - 1.500 rpm	0 - 1.500 rpm	0 - 1.500 rpm
Reversion of rotation direction	—	every 30 s	—	—
Max. length of magnetic bar L x D	30 x 8 mm	30 x 8 mm	30 x 8 mm	30 x 8 mm
Surface plate				
Heating / Surface plate material	polyester foil	polyester foil	glass	glass
Surface plate	115 x 115 mm	Ø = 90 mm	Ø = 110 mm	Ø = 160 mm
General information				
Dimensions W x D x H	115 x 130 x 40 mm	114 x 161 x 12 mm	130 x 135 x 50 mm	180 x 180 x 48 mm
Weight	0,23 kg	0,3 kg	0,48 kg	1 kg
Protection class according to DIN EN 60529	IP 42	IP 65	IP 54	IP 54
Ident-No.	140.2677000	see motifs	see page 20	see page 20