

Magnetomrörare med värme



Magnetic stirrers	RH basic 2	RH basic KT/C safety control	RH digital KT/C safety control	RCT basic safety control
Stirring quantity (H ₂ O)	10 l	15 l	15 l	20 l
Motor rating input / output	15 / 2 W	15 / 2 W	15 / 2 W	12 / 5 W
Speed display	scale (0 - 6)	scale (0 - 6)	digital	scale (1 - 10)
Speed range	100 - 2.000 rpm	100 - 2.000 rpm	100 - 2.000 rpm	0 - 1.100 rpm
Max. length of magnetic bar L x D	40 x 8 mm	40 x 8 mm	40 x 8 mm	80 x 10 mm
Heating function				
Heating output	400 W	500 W	500 W	600 W
Temperature range	RT - 320 °C	RT - 320 °C	RT - 320 °C	RT - 300 °C
Setting accuracy	± 15 K	± 15 K	± 10 K	± 10 K
Adjustable safety circuit	—	100 - 400 °C	100 - 400 °C	50 - 380 °C
Control accuracy with sensor	—	± 3 K	± 3 K	± 1 K
Surface plate				
Heating / Surface plate material	stainless steel	white coated	white coated	silumin
Surface plate	Ø = 125 mm	130 x 130 mm	130 x 130 mm	Ø = 135 mm
General information				
Dimensions W x D x H	168 x 220 x 105 mm	168 x 220 x 105 mm	168 x 220 x 105 mm	160 x 280 x 90 mm
Weight	2,4 kg	2,4 kg	2,4 kg	2,4 kg
Protection class according to DIN EN 60529	IP 21	IP 21	IP 21	IP 42
Interface	—	—	—	—
Recommended accessories (see page 18)	—	No. 2	No. 2	No. 1
Ident-No.	140.3339000	140.3207100	140.3207000	140.3378500