

Multiposition magnetomrörare



Magnetic stirrers	RO 5 / RT 5 power	RO 10 / RT 10 power	RO 15 / RT 15 power
Stirrers	5	10	15
Max. stirring quantity (H ₂ O) per stirrer	0,4 l	0,4 l	0,4 l
Distance between stirring stations	90 mm	90 mm	90 mm
Motor rating input / output	7,2 / 1,8 W	14,4 / 3,6 W	21,6 / 5,4 W
Speed display	scale (1 - 10)	scale (1 - 10)	scale (1 - 10)
Speed range	0 - 1.100 rpm	0 - 1.100 rpm	0 - 1.100 rpm
Deviation for individual stirring stations	5 %	5 %	5 %
Max. length of magnetic bar L x D	30 x 8 mm	30 x 8 mm	30 x 8 mm
Heating function			
Heating output	- / 175 W	- / 375 W	- / 580 W
Temperature range	RT - 120 °C	RT - 120 °C	RT - 120 °C
Setting accuracy	scale (1 - 10)	scale (1 - 10)	scale (1 - 10)
Control accuracy with sensor	± 2 K	± 2 K	± 2 K
Surface plate			
Material	stainless steel / silicone foil	stainless steel / silicone foil	stainless steel / silicone foil
Dimensions	120 x 450 mm	180 x 450 mm	270 x 450 mm
General information			
Dimensions W x D x H	122 x 552 x 65 mm	182 x 552 x 65 mm	272 x 552 x 65 mm
Weight	2,3 - 3 kg	3,2 - 4,2 kg	4,7 - 6 kg
Protection class according to DIN EN 60529	IP 42	IP 42	IP 42
Ident-No.	140.2930200	140.2930400	140.2930600