

506/507



506/507. Midjeförsedda teflonbe- klädda omrörarmagneter

Art.nr	Diam	Längd	Antal/fp
506.20	6	20	5
506.25	6	25	5
506.35	6	35	2
506.45	8	45	2
506.50	8	50	1
506.60	9,5	60	1
506.70	9,5	70	1
507.15019	19	150	1

505



505. Omrörarmagnet, PTFE, slät

Art.nr	Diam	Längd	Antal/fp
505.06	3	6	10
505.08	3	8	10
505.12	4,5	12	10
505.15	4,5	15	10
505.20	6	20	10
505.25	6	25	10
505.30	6	30	10
505.35	6	35	10
505.40	8	40	10
505.45	8	45	10
505.50	8	50	10
505.60	9,5	60	10
505.70	9,5	70	10

502. Omrörarmagnet, PTFE, hantel

Art.nr	Diam	Längd	Antal/fp
502.38	20	38	10
502.56	20	56	10

503. Omrörarmagnet, PTFE, kors

Art.nr	Diam	Längd	Antal/fp
503.06	10	6	10
503.10	20	10	10
503.12	30	12	10

508



508. Omrörarmagnet, PTFE, oval

Art.nr	Diam	Längd	Antal/fp
508.20	10	20	10
508.25	12	25	10
508.30	16	30	10
508.35	16	35	10
508.40	20	40	10
508.50	20	50	10

Utmärkt kemisk, mekanisk och värme
beständighet.
Klarar kontinuerlig arbetstemperaturer
på 250 °C.

Ytterligare varianter på omrörar-
magneter finns

502



503

