

Starkit LM / DV

Innehåller: 3 st pipetter
(1-10ml, 10-100ml, 100-1000ml)
Samt 1 ask pipettspetsar i rack



762. Labmate enkanals
765. Discovery enkanals

Art.nr	Modell	Färgmärkning	Volym (µl)	N (%)	P (%)	Spetsstorlek
762.0002	LM2		Min 0.2	± 12.0	± 6.0	
765.0002	DV2	RÖD	1.0	± 2.7	± 1.3	10 µl
			Max 2.0	± 1.5	± 0.7	
762.0010	LM10		Min 0.5	± 4.0	± 4.0	
765.0010	DV10	RÖD	5.0	± 1.0	± 0.8	10 µl
			Max 10.0	± 0.5	± 0.4	
762.0020	LM20		Min 2	± 4.0	± 3.0	
765.0020	DV20	GUL	10	± 1.0	± 0.7	200 µl
			Max 20	± 0.8	± 0.4	
762.0050	LM50		Min 5	± 2.5	± 2.0	
765.0050	DV50	GUL	25	± 1.0	± 0.6	200 µl
			Max 50	± 0.8	± 0.4	
762.0100	LM100		Min 10	± 1.6	± 0.8	
765.0100	DV100	GUL	50	± 1.0	± 0.4	200 µl
			Max 100	± 0.8	± 0.3	
762.0200	LM200		Min 20	± 1.2	± 0.6	
765.0200	DV200	GUL	100	± 1.0	± 0.3	200 µl
			Max 200	± 0.6	± 0.3	
762.0250	LM250		Min 50	± 1.0	± 0.4	
765.0250	DV250	GRÖN	125	± 0.8	± 0.3	300 µl
			Max 250	± 0.6	± 0.3	
762.1000	LM1000		Min 100	± 1.6	± 0.45	
765.1000	DV1000	BLÅ	500	± 0.7	± 0.25	1000 µl
			Max 1000	± 0.6	± 0.25	
762.5000	LM5000		Min 1000	± 0.6	± 0.3	
765.5000	DV5000	VIT	2500	± 0.6	± 0.3	5000 µl
			Max 5000	± 0.5	± 0.2	
762.X000	LM10000		Min 1000	± 2.5	± 0.6	
765.X100	DV10000	VIT	5000	± 0.8	± 0.3	10000 µl
			Max 10000	± 0.5	± 0.2	

