

560 408-416



560 417-423



562

**560. Normaltratt av propenplast**

Art.nr	Öd	Höjd	Pipud	Pipöd	Piplängd	Förp	Mate.
560.410	31	47	5	6	25	24	PE
560.411	40	65	5	7	34	24	PE
560.412	50	85	9	11	42	24	PP
560.413	75	110	10	13	55	12	PP
560.414	100	155	13	17	80	12	PP
560.415	120	180	14	21	87	12	PP
560.416	150	220	16	24	95	6	PP
560.155	155	215	19	32	115	12	PE
560.417	200	250	20	45	80	6	PE
560.418	250	260	28	48	104	6	PE
560.419	350	440	31	68	175	1	PE
560.422	400	365	42	55	150	1	PP
560.423	450	420	38	55	140	1	PE

911. Pulvertratt av sodaglas

Art.nr	Öd	Höjd	Pip diam
911.080	80	75	24
911.100	100	95	26
911.150	150	135	36
911.200	200	180	38

562. 60°-tratt av propen

Art.nr	Öd	Pipdiam	Piplängd
562.145	25	5	40
562.146	35	6	43
562.147	45	6	45
562.148	65	11	60
562.150	80	12	65
562.152	100	12	80
562.153	120	12	85
562.155	150	16	110
562.156	180	16	145