


Skakapparater med fram och tillbaka rörelse

180

Order order@noaxlab.se



HS 260 control

HS 260 basic

HS 501 digital

Reciprocating shakers	HS 260 basic	HS 260 control	HS 501 digital
Shaking motion	back and forth	back and forth	back and forth
Orbital diameter	20 mm	20 mm	30 mm
Perm. shaking weight with attachment	7,5 kg	7,5 kg	15 kg
Motor rating input / output	45 / 10 W	45 / 10 W	70 / 19 W
Permissible ON-time	100 %	100 %	100 %
Speed range	20 - 300 rpm	10 - 300 rpm	0 - 300 rpm
Speed display	LED-Line	digital	digital
Timer	∞ / 5 - 50 min	∞ / 9 h 59 min	∞ / 0 - 56 min
General information			
Dimensions W x D x H	360 x 420 x 100 mm	360 x 420 x 100 mm	505 x 585 x 120 mm
Weight	8,1 kg	8,8 kg	24 kg
Permissible ambient temperature	5 - 50 °C	5 - 50 °C	5 - 50 °C
Protection class according to DIN EN 60529	IP 21	IP 21	IP 21
Interface	—	RS 232 / analog	—
Recommended accessories (see page 18)	No. 11	No. 11	No. 12
Ident-No.	140.3066600	140.3066700	140.2527000