


Skakapparater med rundgående rörelse

182

Order order@noaxlab.se

MTS 2/4 digital	KS 130 basic	KS 130 control	KS 260 basic	KS 260 control	KS 501 digital
circular	circular		circular	circular	circular
3 mm	4 mm		10 mm	10 mm	30 mm
2 / 4 microtiter plates	2 kg		7,5 kg	15 kg	15 kg
35 / 13,2 W	45 / 10 W		45 / 10 W	70 / 19 W	70 / 19 W
100 %	100 %		100 %	100 %	100 %
0 - 1.100 rpm	80 - 800 rpm		20 - 500 rpm		0 - 300 rpm
scale	LED-Line	digital	LED-Line	digital	digital
∞ / 1 - 99 min	∞ / 5 - 50 min	∞ / 0 - 9 h 59 min	∞ / 5 - 50 min	∞ / 0 - 9 h 59 min	∞ / 0 - 56 min
185 x 320 x 105 mm	270 x 316 x 98 mm		360 x 420 x 98 mm		505 x 585 x 120 mm
2,7 kg	7,7 kg	8,8 kg	8,1 kg	8,8 kg	24 kg
5 - 40 °C	5 - 50 °C		5 - 50 °C		5 - 50 °C
IP 21	IP 21		IP 21		IP 21
—	RS 232 / analog		RS 232 / analog		—
incl.	No. 10		No. 11		No. 12
140.3208000	140.2980000	140.2980100	140.2980200	140.2980300	140.2526400