

486

**486. Burk och lock av solid PTFE**

Art.nr	Volym	Diam	Höjd	Hals
486.0100	100	52,56	89,82	35,24
486.0250	250	66,74	118,40	42,00
486.0500	500	80,00	147,15	51,15
486.1000	1000	102,00	184,21	57,00

Utmärkt kemisk, mekanisk och värme beständighet.
Klarar kontinuerlig arbetstemperatur på 250 °C


515. Teflonproppar, NS

Art.nr	NS	Förp
515.10	10/19	2
515.14	14/23	2
515.19	19/26	2
515.24	24/29	1
515.29	29/32	1
515.34	24/35	1

515



116

**116. Skopa, PTFE**

Art.nr	Volym	Diam	Höjd
116.0010	10	25	34
116.0100	100	55	65
116.0250	250	70	95
116.0500	500	80	125
116.1000	1000	105	155

95

Order tel 08-581 700 55 fax 08-581 700 65